ANTIGENOS FEBRILES: PORTA

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PRO4_REG9_AFPORTA

FEBRILE ANTIGENS: SLIDE

QUALITATIVE AND SEMIQUANTITATIVE DETERMINATION

For "in vitro" diagnostic

Principle

Febrile antigens are bacterial suspensions specially prepared for use in slide agglutination procedure for detecting specific agglutinins implicated in certain febrile diseases such as Brucellosis and Salmonellosis.

Reagents

S. Typhi H (FlagellarAg., d-H).)	5 ml	Ref. 99 01 88
S. Typhi O (Somatic Ag., 9, 12-O)	5 ml	Ref. 99 81 01
S. Paratyphi A H (Flagellar Ag., a-H)	5 ml	Ref. 99 99 17
S. Paratyphi B H (Flagellar Ag, b-H.)	5 ml	Ref. 99 08 70
Brucella Abortus	5 ml	Ref. 99 90 01
Brucella Melitensis	5 ml	Ref. 99 82 79
Proteus OX 19	5 ml	Ref. 99 61 15

Controls

Salmonella Positive Control	1 ml	Ref. 99 91 29
Brucella Positive Control	1 ml	Ref. 99 46 11
Proteus Positive Control	1 ml	Ref. 99 91 31
Negative Control	1 ml	Ref. 99 39 10

Storage and stability

When stored at 2 - 8° C, the reagents will remain stable until the expiration date stated on the label.

Caution

Somatic antigens contain 0,5% phenol and flagellar antigens contain 0,5% formalin as preservatives.; they should be handled with care
 Waste products must be handled as per local regulations

Sample

Serum. Can be stored up to 1 week at 2 - 8° C.

A. Slide agglutination. Screening test

Procedure
1.- Bring sample, reagent and controls to room temperature.
2.- Softly homogenize the bacterial suspension including the volume in the dropper.
3.- Dispense 20 µl of serum and each control on different zones of the slide. Add one drop of the bacterial suspension to each zone and mix with the aid of a disposable stirrer.
4.- Tilt during 1 min. and observe with the aid of a proper light, for the presence or absence of agglutination.
Results
POSITIVO: Agglutination. Equivalent to a titre of 1:80 (aprox) by the Tube technique.
NEGATIVO: No agglutination. Titre lower than 1:80.

B. Slide agglutination. Semicuantitative determination

Procedure
1.- Bring sample, reagent and controls to room temperature.
2.- Softly homogenize the bacterial suspension including the volume in the dropper.
3.- Dispense 80, 40, 20, 10 and 5 µl of serum and each control on different zones of the slide. Add one drop of the bacterial suspension to each zone and mix with the aid of a disposable stirrer.
4.- Tilt during 1 min and observe with the aid of a proper light, for the presence or absence of agglutination.
Results
The serum volumes dispensed above correspond approximately to titres of 1:20, 1:40, 1:80, 1:160 and 1:320 by the Tube technique.
POSITIVO: Agglutination.
NEGATIVO: No agglutination.

Interpretation of results

Titres equal or higher than 1:80 have probable diagnostic significance. However, the precise diagnosis of the disease depends not only on the analytical results but also on a whole clinical evaluation. Titre rises considerably between the acute phase and the convalescence.

Performance characteristics

The antigenicity of all febrile antigens have been standardized with control antiserum of established titer

Interferences: There is no interference by RF at titers lower than 300 IU/ml. Hemoglobin (≤10 g/L/L); Lipids (≤10 g/L) or Bilirubin (≤20 mg/dl) do not interfere the test

Traces of detergent in the labware may give falsely positive results, as well as the use of an old or contaminated saline solution.

QUIMICA CLINICA APLICADA S.A.
ISO 9001 / ISO 13485 Certified Company
 A7 Km 1081 – P.O. Box 20 - E43870 AMPOSTA / SPAIN
 Tel. ++ 34 (977) 70. 62. 30 Fax ++ 34 (977) 70. 30. 40
Revision: October 2011

Negative results can be obtained with patients at the early stages of the disease or with those under antibiotic treatment.

In some cases, prozone phenomena can be found; thus, when suspecting a high titre, it is advisable to dilute the sample to 1:20 with physiological saline (NaCl 0.9 %) and perform the test once again.

Quality control

It is recommended that positive and negative controls be run with the samples to be processed to check the right performance of the test.

C.- Tube agglutination. Quantitative determination

It is recommended to confirm each positive reaction by the **Tube technique**.

The test can be done with the same reagents used for slide test as it is stated on the following method

Procedure	Reading
For each antigen prepare, as indicated in the following scheme 7 Kahn tubes and place them in a rack.	Using a proper lamp and reading against a black background, record the degree of agglutination as indicated:
Tube Nº. 1 2 3 4 5 6 7	++++ 100% of the microorganisms are agglutinated. Supernatant is clear.
Saline ml 1,9 1,0 1,0 1,0 1,0 1,0 1,0	
Serum ml 0,1 1,0 1,0 1,0 1,0 1,0 1,0	+++ Approximately 75% of the microorganisms are agglutinated. Supernatant is clear
Antigen µl 50 50 50 50 50 50 50	++ Aproximately 50% of the microorganisms are agglutinated. Supernatant is moderately cloudy.
Título 1:20 1:40 1:80 1:160 1:320 1:640 1:1280	+ Approximately 25% of the microorganisms are agglutinated. Supernatant is cloudy.
1.- Pipette 1.9 ml of saline into the first tube of each row and 1.0 ml into the remaining tubes.	- No agglutination. Microorganisms remain as a cloudy suspension.
2.- Add 0.1 ml of the serum to be assayed into the first tube, mix well, and transfer 1.0 ml to the second tube.	
3.- Mix well and transfer again 1.0 ml from tube Nr.2 to tube Nr.3 and so on. Discard 1.0 ml from tube Nr.7 after mixing thoroughly.	
4.- Add 50 µl of the desired antigen to each of the 7 tubes shake the rack and incubate at 37° C for 24/48 h.	

Materials required but not provided

For A and B methods
 Glass or plastic slides
 Disposable stirrers

For C method:
 Khan tubes
 Thermostated water bath or similar at 37°C

For all methods
 Pipettes
 Saline solution (NaCl 0.95)

NOTE

The tone and colour intensity of the antigen suspension can vary from lot to lot. This variation will not affect the final results.

Quantitative determination

QCA has also available the following reagents, diluted and ready for use:

S. Typhi H (Ag flagelar, d-H)	90 ml	Ref. 99 63 92
S. Typhi O (Ag somático, 9,12-O)	90 ml	Ref. 99 90 19
S. Paratyphi A H (Ag flagelar,a -H)	90 ml	Ref. 99 17 93
S. Paratyphi B H (Ag Flagelar, b-H)	90 ml	Ref. 99 55 76
Brucella	90 ml	Ref. 99 01 42

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